

ABRF-LMRG - 3D Microsphere Standard Sample Study #3

Imaging Protocol FV1000 – 3D Point Spread Function (Slide1)

- 1) Turn on and warm up the lasers for at least one hour. Take the DIC elements OUT of the lightpath.
- 2) Put a 1.0 NA or higher objective lens in place. If necessary, put immersion media on the lens. (Recommend using the Olympus 60x oil objective for one of your tests, if available).
- 3) Load Alexa 488 in the Dye List to facilitate setting up the lightpath quickly. The **DM 405/488** nm Main Beam Splitter will be selected. Using the VBF (Variable Barrier Filter), adjust the emission bandwidth of Channel 1 to **~495-540 nm**.
- 4) In the **Image Acquisition Control** window, set the sequential mode to “**Line**”, **choose a Kalman of 4**. Set up ‘**Harddisk recording**’. In the **Acquisition Setting** menu, moving from the top to bottom of the window: Choose **one-way imaging** (left to right), choose a pixel dwell of **2 us/pixel**, choose a **Frame Size** of **512x512 pixels**, choose a **Zoom factor** to achieve a final pixel size of **0.07 μ m** (see ‘i’ button under the Image Acquisition Control window for pixel size information). Choose a z step size of **0.1** (for PSF tests).
- 5) Set the **Pinhole** to 1 Airy unit (AU) which is chosen when ‘Auto’ is selected under the Confocal Aperture.
- 6) Set each detection channel photomultiplier tube (PMT) - indicated by **HV** under the Image Acquisition control menu - to 700-800 units. Use **Ctrl H** (Hi-Lo) to check the saturation levels. An average intensity of ~3500 units for each colour channel is ideal but it is important to ensure that no pixels are saturated with these gain settings. **Note:** there may be microsphere aggregates larger than 0.1 μ m in your sample; do not include these when setting the intensity thresholding.
- 7) The **Digital Offset** should be set to zero – so that no pixels read zero intensity units. Double check there are no zero intensity (blue) pixels using the CTRL H (Hi-Lo) feature.
- 8) Set the **Digital Gain** to 1.
- 9) Suggested laser intensity is: Green–25 mW Argon Ion-488 nm laser = 1%. Recheck the image and ensure that intensity values of ~3500 grey levels are maintained with these specific laser settings and that no pixels measure saturated intensities.
- 10) Press the **XY button** (no depth) to take an image of the microspheres in the plane of focus.
- 11) Perform a final verification of the image acquisition settings using **CTRL H** for the plane of focus. Zero intensity pixels will be displayed as blue and saturated pixels as red. If there are blue pixels, increase the **Digital Offset**. If there are red pixels reduce the **laser power (not the HV setting!)**.
- 12) Choose the **Z-Stack** option and set up the Z-axis scanning with a total z-stack size of 100 μ m and a step size of 0.1 μ m between images. Use the **Live** scanning mode and set the **Z-Stack** using the **Start/End** menu. Focus below the sample and mark the **first plane** after the sample goes out of focus as 0. Set the **last plane** at 100 μ m above the **first plane**. The image should be automatically saved to harddisk recording since it will be a large stack.
- 13) Press the **XYZ** button to perform the **Z-Stack** acquisition.
- 14) Save the image stacks as .oib files and also as .tif files. Name the files as follows:

IMPORTANT:

“**Lastname_Firstname_Microscope_Platform_MagnificationX_NA_ImmersionMedia_PSF**”
e.g. **Brown_Claire_Zeiss710_63X_NA1.4_Oil_PSF**)

- 15) Please follow this link to submit a sample information form for each dataset you submit
<https://tinyurl.com/LMRG-vial1> (takes a few seconds to load correct link if you click on it).

ABRF-LMRG - 3D Microsphere Standard Sample Study #3
Imaging Protocol FV1000 – 3D Point Spread Function (Slide1)