

# ABRF-LMRG - 3D Microsphere Standard Sample Study #3: SNR

## Imaging Protocol: Nikon A1R – 3D Signal-to-Noise (Slide 2)

### This protocol is for the SNR acquisition only

- 1) Turn on and let your lasers warm up for at least 1 hour.
- 2) Put a **0.5 NA** or higher objective lens in place. If necessary, put immersion media on the lens.
- 3) From the **OC** panel, select the appropriate Optical Configuration(s) that will allow you to image:
  - a. Green Fluorescent Microspheres: 488nm ex, 505nm em
  - b. Red Fluorescent Microspheres: 555nm ex, 575nm em
- 4) Within the **A1plus Compact GUI** panel, select the following imaging conditions:
  - a. Galvano
  - b. Unidirectional scan
  - c. Pixel Dwell = 2.2  $\mu$ s
  - d. Size = 512
  - e. Pinhole = 1.2 AU (calculated for 488nm)
  - f. Line Average = 4x
  - g. Channel Series
- 5) For the each laser, set the illumination power to **1%**, and the **HV gain (HV(G))** for the PMT to **30** units. Finally, set the **Offset** to **0**.

**Note:** If you do not have GaAsP detectors, you may need to increase your **HV(G)** to 90. If you are unsure what detectors you do have, start at 30 and work your way up during a live scan.
- 6) Start a live scan and find a viable imaging region. An ideal region will have many beads in the field of view, but separate enough to generate distinct beads. Bring the beads into focus.
- 7) Select the **Pixel Saturation Indication** icon and check for saturated pixels. Set the Pixel Saturation Indication color to one that can be clearly observed on both laser channels (Complimentary Color tends to work well).
- 8) Adjust your **laser power** and your **HV(G)** to avoid saturation while generating a peak pixel intensity value of approximately **3500 counts**. Check your settings by scrolling through multiple z planes as the beads vary in intensity.
- 9) In the **A1plus Scan Area** tab, select a square scan area (first **icon** on the top left, in frame scan mode). Choose a **Pixel size** to **0.2  $\mu$ m** per pixel (or as close as possible).
- 10) Set up the acquisition of Z stacks within the **ND Acquisition** window (or press Ctrl+Alt+shift+Z). Lower your objective to a focal plane just below the initial layer of microspheres. Choose the **Asymmetric** option within the **Z stack** tab within the **ND Acquisition** Window. Set the current focal plane to home by selecting the **Home** icon. Set **Below** as 0 and **Above** as +100  $\mu$ m. Set the step size to **1.0  $\mu$ m**. There should be 101 z steps.
- 11) Press the **Run Now** button to perform the acquisition.
- 12) Please follow this link to submit a sample information form for each dataset you submit:  
<https://tinyurl.com/LMRG-vial1>
- 13) Save the files as **.nd2** and also as **.tif**. Name the files as follows:  
“**Lastname\_Firstname\_Microscope\_Platform\_MagnificationX\_NAY\_ImmersionMedia**”.  
For example: **Brown\_Claire\_NikonA1R\_100x\_NA1.45\_Oil**