

ABRF-LMRG - 3D Microsphere Standard Sample Study #3: SNR

Imaging Protocol: Nikon C2 – 3D Signal-to-Noise (Slide 2)

This protocol is for the SNR acquisition only

- 1) Turn on and let your lasers warm up for at least 1 hour.
- 2) Put a **0.5 NA** or higher objective lens in place. If necessary, put immersion media on the lens.
- 3) From the **Dye and Spectral Settings** panel, select the appropriate Optical Configuration(s) that will allow you to image:
 - a. Green Fluorescent Microspheres: 488nm ex, 505nm em
 - b. Red Fluorescent Microspheres: 555nm ex, 575nm em
- 4) Within the **A1plus Compact GUI** panel, select the following imaging conditions:
 - b. Unidirectional scan
 - c. Pixel Dwell = 1.9 μ s
 - d. Size = 512
 - e. Pinhole = 1.2 AU (calculated for 488nm)
 - f. Line Average = 4x
 - g. Channel Series
- 5) For the each laser, set the illumination power to **1%**, and the **HV gain** for the PMT to **90** units. Finally, set the **Offset** to **0**.
- 6) Start a live scan and find a viable imaging region. An ideal region will have many beads in the field of view, but separate enough to generate distinct beads. Bring the beads into focus.
- 7) Select the **Pixel Saturation Indication** icon and check for saturated pixels. Set the Pixel Saturation Indication color to one that can be clearly observed on both laser channels (Complimentary Color tends to work well).
- 8) Adjust your **laser power** and your **HV** to avoid saturation while generating a peak pixel intensity value of approximately **3500 counts**. Check your settings by scrolling through multiple z planes as the beads vary in intensity.
- 9) In the **C2plus Scan Area** tab, select a square scan area (first **icon** on the top left, in frame scan mode). Choose a **Pixel size** to **0.2 μ m** per pixel (or as close as possible).
- 10) Set up the acquisition of Z stacks within the **ND Acquisition** window. Lower your objective to a focal plane just below the initial layer of microspheres. Choose the **Asymmetric** option within the **Z stack** tab within the **ND Acquisition** Window. Set the current focal plane to home by selecting the **Home** icon. Set **Below** as 0 and **Above** as +100 μ m. Set the step size to **1.0 μ m**. There should be 101 z steps.
- 11) Press the **Run Now** button to perform the acquisition.
- 12) Please follow this link to submit a sample information form for each dataset you submit:
<https://tinyurl.com/LMRG-vial1>
- 13) Save the files as **.nd2** and also as **.tif**. Name the files as follows:
“**Lastname_Firstname_Microscope_Platform_MagnificationX_NAY_ImmersionMedia**”.
For example: **Brown_Claire_NikonA1R_100x_NA1.45_Oil**