

ABRF-LMRG - 3D Microsphere Standard Sample Study #3

Imaging Protocol Zeiss 710 – 3D Signal-to-Noise (Slide 2)

- 1) Turn on and warm up the lasers for at least one hour.
- 2) Put a 0.5 NA or higher objective lens in place. If necessary, put immersion media on the lens.
- 3) Under the **Light Path** window set up the confocal for imaging green and orange. For example under **Visible Light**, set up the light path using the 488 nm laser and either a 543 nm or 561 nm laser. Use the **488/543 nm Main Beam Splitter (MBS 488/543)** or the **488/561 nm Main Beam Splitter (MBS 488/561)** to reflect the two laser excitation lines. Set up two channels with emission bandwidths of **~495-540 nm** and **~560-650 nm**.
- 4) In the **Acquisition Mode** window set the scan mode to **Frame**. Choose a **Frame Size** of 512x512 pixels, **Line Step** of **1**, **Scan Speed** of **9** (pixel dwell time of **~1.27 µs/pixel**), **Mean Line Averaging** of **4** to reduce pixel noise, **Bit Depth** of **12 Bit** and adjust the **Zoom** factor to achieve a **Pixel Size** between **0.20 and 0.30 µm**. Under **Direction** set for unidirectional scanning (indicated by right pointing arrow).
- 5) Set the **Pinhole** for each channel to 1 Airy unit (AU).
- 6) Set each detection channel photomultiplier tube (PMT) or **Gain (Master)** to 700-800 units. An average intensity of ~3500 units is ideal for each colour channel but it is important to ensure that no pixels are saturated with these gain settings.
- 7) The **Digital Offset** must be set above zero (~12) so that no pixels read zero intensity units.
- 8) Set the **Digital Gain** to 1.
- 9) Suggested laser intensities are: Green–25 mW Argon Ion-488 nm laser = >1% and Orange–20 mW-561 nm laser = >1%. Recheck the images for each channel and ensure that intensity values of ~3500 grey levels or maintained with the specific laser settings and that no pixels measure saturated intensities.
- 10) Press the **Snap** button to take an image of the microspheres in the plane of focus.
- 11) Perform a final verification of the image acquisition settings for each channel using the **Range Indicator** Look-Up Table (LUT). Click on the colour bar for each track in the image to change the colour. Under the **Dimensions** tab at the bottom of the image select the palette and choose the **Range Indicator** LUT. Zero intensity pixels will be displayed as blue and saturated pixels as red. If there are blue pixels increase the **Digital Offset**, if there are red pixels reduce the **laser power**.
- 12) Choose the **Z-Stack** option and set up the Z-axis scanning with a total z-stack size of 100 µm and a step size of 1.0 µm between images. Use the **Live** scanning mode and set the **Z-Stack** using the **First/Last** window. Focus below the sample and mark the **first plane** after the sample goes out of focus. Set the **last plane** at 100 µm above the **first plane**.
- 13) Press the **Start Experiment** button to perform the **Z-Stack** acquisition.
- 14) Save the image stacks as .ism files and also as .tif files. Name the files as follows:

IMPOTANT:

“**Lastname_Firstname_Microscope_Platform_MagnificationX_NAY_ImmersionMedia_SNR**”
e.g. **Brown_Claire_Zeiss710_63X_NA1.4_Oil_SNR**)